



New Q-Band Artificial Intelligence - Based Chromosome Analysis and Karyotyping a Pilot Evaluation

Liz Pardo, BSc, MSc, PhD(c)¹, Ivan Acosta, BSc¹, Claudia Acosta, BSc¹, Carolina Villamizar, BSc¹,
Cristina Steele, CG (ASCP), MBA², Einat Zindany, BSc², Yael Glickman, PhD, MBA²

¹ Cytogenomics Laboratory - Genomics Sciences Unit, COLCAN-ANNAR; ² Applied Spectral Imaging, Carlsbad, CA, USA

Introduction

- The **Q-banding karyotyping method** is broadly used to identify individual chromosomes and their structural anomalies. In this technique, Quinacrine stain causes the AT-rich regions of the chromosomes to fluoresce more brightly than the GC-rich regions, therefore creating Q-bands.
- Recent advancements in artificial intelligence (AI) have been integrated into **G-banded karyotyping platforms**, significantly reducing the number of manual interventions required to correct a metaphase when compared to conventional digital imaging systems that rely on traditional image processing technologies.
- Multicenter studies performed on G-banded specimens have reported **91% correct chromosome segmentation** and **97% correct chromosome placement** when using new AI capabilities. (1-5)
- The present pilot evaluation reports preliminary assessment of a new **AI-based Q-banded karyotyping** application.




Figure 1 Representative examples of automatically segmented metaphases (top) and automatically placed chromosomes in karyograms following manual segmentation adjustments (bottom)
Left: bone marrow. Right: peripheral blood. Errors are highlighted).

Methods

- Standard cytogenetic cell cultures were processed and spread onto microscope slides for staining. **Q-banding slides were scanned and analyzed** using the **HiBand system** (Applied Spectral Imaging).
- Metaphases, automatically identified by the system at 10X magnification, were captured at 100X. A newly developed AI-based computer-aided karyotyping application was used to automatically analyze selected cells.
- System-suggested aided karyograms were generated based on a predefined range of detected objects (35-49).
- Aided karyograms were compared with the karyotype table prepared by certified cytogenetics professionals to assess the number of correctly segmented and classified chromosomes.
- The number of manual adjustments required to correct the automatically suggested karyograms were then evaluated.

Results

- Ten Q-band peripheral blood and ten Q-band bone marrow metaphases** with an average of 9 ± 6 overlaps and/or touching chromosomes were included in this pilot evaluation.
- Aided karyograms were generated for **85% of the metaphases**, representing a total of 782 chromosomes.
- The system **automatically segmented 94%** of the blood chromosomes (390 out of 413) and **79%** of bone marrow chromosomes (292 out of 369).
- Following segmentation correction, 408 blood chromosomes (**99%**) and 351 bone marrow chromosomes (**95%**) were **correctly placed** in the karyotype table.
- For both cell types, 98% of the correctly placed chromosomes were also correctly oriented.
- In average, **3±2 manual adjustments** were necessary to correct both segmentation and classification of peripheral blood metaphases, and **10±6 manual adjustments** were necessary to correct bone marrow metaphases.




Figure 2 Representation of the Q-banded karyotyping workflow
* Following manual adjustments to correct automatic segmentation

Discussion and Conclusions

- This pilot evaluation is the first report of a new AI-powered software for the analysis of Q-banded stained slides.
- The AI software **correctly segmented 87%** of peripheral blood or bone marrow chromosomes in metaphases containing an average of 9 overlaps and/or touching chromosomes.
- In addition, **97% of the chromosomes were correctly placed** in the karyotype table following segmentation correction.
- These results are **aligned with accuracy reported for G-banded specimens**.
- This pilot evaluation should be expanded to further confirm these preliminary results.

Disclosures

LP, IA, CA and CW have no disclosure. CS, EZ, YG are employees of Applied Spectral Imaging.

References

- Taghiyev A et al. Genetics in Medicine Open 2024, 2(1):101575
- Burnside R et al. Genetics in Medicine Open 2025, 3(2):103053
- Burnside R et al. Genetics in Medicine Open 2025, 3(1):101918
- Steele C et al. Genetics in Medicine Open 2025, 3(2):103033
- Rosenblum L et al. Genes 2025, 16(6):285

